

Helping people prevent, prepare for,
and respond to emergencies.



2009-2010 Training Products Catalog



You are a valued customer of the American Red Cross, and we are pleased to provide you with this catalog featuring quality products to support your training needs. Our philosophy is simple: American Red Cross' Quality Training and Laerdal's Quality Products make for Quality Care.

We appreciate your commitment in promoting lifesaving skills. We believe that by working together at all levels we can assure that you are able to provide those around you with the skills necessary to save a life.

The American Red Cross constantly strives to respond to the health and safety concerns of Americans at home, in school and in the workplace. Working together with you, we can all meet the challenges of today and be prepared for tomorrow. Thank you for all that you do!

“Laerdal and the American Red Cross have provided me with so much useful information and good products.

I don't have to shop around.”

- Rose Aldrich, Authorized Provider,
MD Medical Center

Table of Contents 1-2

American Red Cross Health and Safety Courses..... 3-4

First Aid/CPR/AED Training 5-12


CPR Breathing Barriers 13-16

Aquatic Training 17

Professional Rescuer Training and Equipment 18-33

Ordering Information and Customer Support 34-35





The American Red Cross offers a variety of health and safety training opportunities. If you are interested in any of the following area please contact your local Red Cross Chapter.

Course Options for the Workplace

- First Aid/CPR/AED for the Workplace*
 - Adult CPR/AED*
 - Adult CPR*
 - Adult AED*
- Infant and Child CPR/AED*
- First Aid and Preparedness
- Bloodborne Pathogens:
 - Preventing Disease Transmission
- Administering Emergency Oxygen
- Injury Control Modules:
 - Ergonomics
 - Slips, Trips, and Falls
 - Back Injury Prevention
 - Workplace Violence Awareness
 - Managing Stress
 - Your Heart Matters/Act in Time

Course Options for the Community

- First Aid/CPR/AED for Schools and the Community*
 - First Aid*
 - Adult CPR/AED*
 - Adult CPR*
 - Adult AED*
 - Infant and Child CPR/AED*
- Swimming and Water Safety
 - Learn-to-Swim
 - Longfellow's Whales Takes
 - Guard Start: Lifeguarding Tomorrow
- Family Caregiving
- Nurse Assistant Training
- Babysitter's Training
- Sport Safety Training
- Basic Aid Training (B.A.T.)
- First Aid for Children Today (F.A.C.T.)
- Pet First Aid

Course Options for Professional Rescuers

- Emergency Response
- CPR/AED for the Professional Rescuer
- Lifeguarding
- Bloodborne Pathogens:
 - Preventing Disease Transmission
- Administering Emergency Oxygen

For more information on the Red Cross courses listed here, visit RedCross.org

The American Red Cross always needs enthusiastic and committed individuals to teach its programs. For information on how to become an American Red Cross instructor or instructor's aide, contact your local chapter or visit RedCross.org.

Course availability varies between chapters. Contact your local chapter for a full course listing in your area.

*Available in Spanish



Laerdal® Little Anne® Manikin

Realistic and inexpensive quality adult CPR training

The Little Anne manikin is a realistic, inexpensive, lightweight manikin that offers all the essential features for CPR training. Designed with both the user and instructor in mind.

- Airway can be opened using head-tilt/chin-lift or jaw-thrust technique
- A “clicker” confirms correct compression depth
- Compression spring provides consistent resistance
- Economical, disposable airways equipped with non-rebreathing valve
- Removable, reusable faces



All Little Anne Carry Cases are Red Cross Branded

- 02 00 21 Little Anne Manikin
- 02 00 26 Little Anne Manikin (Dark Skin)
- 02 00 23 01 Little Anne 4-Pack
- 02 00 28 01 Little Anne 4-Pack (Dark Skin)

All Little Anne Manikins Include:

Single Pack: 1 Manikin, Soft Pack/Training Mat, 2 Resusci Manikin Faces, 2 Little Anne Airways, 6 Manikin Wipes, and Directions for Use.

*Four Pack: 4 Manikins, 4 Soft Pack Carry Case/Training Mat, 8 Resusci Manikin Faces, 8 Little Anne Airways, 24 Manikin Wipes, and Directions for Use. All in a 4-Pack Carry Case.**

**4-Pack Carry Case is equipped with wheels for easy portability.*



For CPR breathing barriers, see pages 13-16.



Baby Anne® Infant CPR Trainer

Realistic and inexpensive quality infant CPR training

- Lightweight, portable four pack weighs under 12 pounds
- Affordable price for increased hands-on classroom practice
- Realistic chest compression and rise
- Inexpensive and easy-to-maintain airway system (one airway per student)
- Unique foreign-body airway obstruction feature
- Head tilt/chin lift is required to open the airway



Baby Anne features a realistic chest compression and rise.



Release of foreign-body airway obstruction can be practiced.



All Baby Anne Carry Cases are Red Cross Branded

All Baby Anne Manikins Include:

Single Pack: 1 Manikin, Soft Pack Carry Case, 6 Airways, 10 Foreign-Body Practice Objects, and Directions for Use.

Four Pack: 4 Manikins, Soft Pack Carry Case, 24 Airways, 40 Foreign-Body Practice Objects, and Directions for Use.

- 05 00 01 Baby Anne Manikin
- 05 00 03 Baby Anne Manikin (Dark Skin)
- 05 00 11 Baby Anne Manikins 4-Pack
- 05 00 13 Baby Anne Manikins 4-Pack (Dark Skin)

Baby Anne can also be used in your American Red Cross Babysitter's Training classes!

Little Junior™ Child CPR Trainer

Realistic and inexpensive quality pediatric CPR training

The Little Anne manikin is a realistic, inexpensive, lightweight adult CPR trainer that offers all the essential features. Designed with both the user and instructor in mind.

- Airway can be opened using head-tilt/chin-lift or jaw-thrust technique
- Compression “clicker” provides audible confirmation of adequate compression depth
- Economical, disposable airways equipped with non-rebreathing valve for quick and easy cleanup
- Removable, reusable faces

- 18 00 21 01 Little Junior Manikin
- 18 00 26 01 Little Junior Manikin (Dark Skin)
- 18 00 23 01 Little Junior 4-Pack
- 18 00 28 01 Little Junior 4-Pack (Dark Skin)

All Little Junior Manikins Include:

Single Pack: 1 Manikin, 1 Soft Pack Carry Case/Training Mat, 2 Faces, 2 Airways, 6 Manikin Wipes, and Directions for Use.

Four Pack: 4 Manikins, Soft Pack Carry Case/Training Mat, 8 Junior Faces, 8 Junior Airways, 24 Manikin Wipes, and Directions for Use.



All Little Junior Carry Cases are Red Cross Branded


Little Family Pack

Now Laerdal brings you realistic, affordable, quality CPR training in one “Little” happy family package – the new Little Family Pack. Designed to reflect the realistic anatomical differences among an adult, a child, and an infant.

- Complete with one Little Anne®, one Little Junior, and one Baby Anne®, all in their individual cases.

- 020081 Little Manikin Family Pack,ARC
- 020086 Little Manikin Family Pack Dark Skin,ARC

Little Family Pack Carry Cases are not Red Cross Branded



Costs less than buying all three manikins separately. It's like getting Baby Anne FREE!

Airways and Faces

- Little Anne®
- 02 03 00 Little Anne Airways (Pkg. 24)
 - 02 03 01 Little Anne Airways (Pkg. 96)
 - 31 02 10 Manikin Faces (Pkg. 6)
 - 31 02 15 Manikin Faces, Dark Skin (Pkg. 6)

- Baby Anne™
- 05 01 00 Airways (Pkg. 24)
 - 05 02 00 Baby Anne Faces (Pkg. 6)
 - 05 02 02 Baby Anne Faces, Dark Skin (Pkg. 6)

- Little Junior™
- 18 30 10 Junior Faces (Pkg. 6)
 - 18 30 15 Junior Faces, Dark Skin (Pkg. 6)
 - 18 32 10 Junior Airways (Pkg. 25)
 - 18 32 11 Junior Airways (Pkg. 100)

Baby Anne Accessories

- 05 03 00 Foreign-Body Practice Objects (Pkg. 10)
- 05 04 00 Jaw
- 05 05 00 Mask Connector
- 05 07 00 Overalls with Breastplate

Miscellaneous

- 15 24 00 Manikin Wipes (Pkg. 50)
- 15 24 01 Manikin Wipes (Pkg. 1,200)
- 15 12 01 03 Manikin Face Shields (6 Rolls of 36 Each)



Red Cross Branded Carrying Case

- Little Anne
- 02 08 00 Soft Pack/Training Mat
 - 02 08 10 Soft Pack/4-Pack

- Little Junior
- 18 38 21 Soft Pack/Training Mat
 - 18 38 31 Soft Pack/4-Pack

- Baby Anne
- 05 06 01 Soft Pack
 - 05 06 11 Soft Pack/4-Pack

The American Red Cross Automated External Defibrillation (AED) Program

Sudden Cardiac Arrest (SCA) is a leading cause of death in the United States, claiming approximately 300,000 lives every year. Up to 50,000 of these deaths could have been prevented if someone on scene had initiated the Cardiac Chain of Survival and an automated external defibrillator (AED) had been available for immediate use at the time of the emergency.

The American Red Cross has a vision of one person in every household being trained in First Aid and CPR lifesaving skills, and all Americans being within four minutes of an AED and someone trained to use it in the event of Sudden Cardiac Arrest.

To meet this vision, the Red Cross is helping communities plan and implement access to defibrillation programs by providing AED education and training, supporting state and federal AED legislation, as well as urging communities to obtain and place AEDs in public places.

The American Red Cross Can Facilitate Your Purchase of AEDs

To help meet its mission of saving lives, the Red Cross has formed relationships with 4 AED manufacturers (Cardiac Science, Philips Health Care Group, Physio-Control, a division of Medtronic, and ZOLL Medical Corporation) to facilitate the purchase of AED devices for customers. If your company is interested in developing a workplace AED training program or obtaining an AED, contact your local American Red Cross chapter or visit RedCross.org. A portion of the purchase price goes back to your local chapter to support services and activities in your community.



AED Little Anne® Training System

The AED Trainer 2 with remote control, AED Little Anne manikin, and special Laerdal LINK training pads work together to provide the student with CPR-D feedback

Special Features

- Chest skin area illuminates when pads are positioned properly. With no visible connectors on the manikin, the student must utilize anatomical landmarks for proper pad placement. Reinforces proper AED training skills
- Feedback is provided via voice prompts from the AED Trainer 2 and LEDs on the manikin
- Pads must be properly placed on the chest skin for the AED Trainer 2 to proceed with its scenarios
- 10 pre-programmed training scenarios included
- Add the Programming Kit to design your own customized training scenarios

02 00 62 Little Anne AED Training System

02 00 52 AED Little Anne (without AED Trainer 2)

System includes:

AED Little Anne, Laerdal AED Trainer 2 with Soft Pack, Carry Case for Little Anne, Remote Control, and Laerdal LINK Training Pads.



Also available as separate products, see pages 11 and 19.



02 50 31 American Red Cross Carry Case

02 50 21 AED Little Anne with American Red Cross Jacket

94 50 90 Training Pads, Laerdal LINK Training System

94 50 32 American Red Cross Soft Case

94 00 52 01 Laerdal AED Trainer 2

94 50 50 Remote Control for Laerdal AED Trainer 2

Laerdal® AED Trainer 2

The Laerdal AED Trainer 2 provides a realistic and affordable AED training solution by simulating the HeartStart® FR2 AED in appearance, operation, and voice prompts

- 10 pre-programmed scenarios
- Standard training pads for use with any CPR manikin
- Status indicator window can be changed manually by the instructor
- Includes high-quality, zippered soft pack carry case
- New battery pack simulates the FR2 battery (6 C-Cell batteries included)



Laerdal AED Trainer 2 Includes:
Trainer, Soft Pack, 6 Batteries, Standard Adult Training Pads.

Battery Pack	Remote Control	Programming Kit	Soft Pack

- 94 00 52 01 Laerdal AED Trainer 2 (includes Trainer, American Red Cross Soft Pack Carry Case, 6 C-Cell Batteries, 1 Set of Standard Training Pads)
- 94 50 50 Remote Control
- 94 50 55 Programming Kit

HeartStart® AED Trainer

Designed for realistic and specific AED training

- Eight CPR/AED pre-configured scenarios
- CPR coaching feature
- Offers step-by-step instructions for performing CPR, including compression-to-ventilation ratios and the length of each ventilation
- Simulates operation of the HeartStart OnSite or Home Defibrillators only



- M5085A-ABA HeartStart AED Trainer
- M5094A Pediatric Training Pads

HeartStart AED Trainer Includes:
Trainer, Soft Pack Carry Case, Adult Training Pads Cartridge (1), External Manikin Adaptor (1), Directions for Use and a Quick Reference Guide.

AED Training Accessories

AED Little Anne® Training System

- 02 50 31 American Red Cross AED Little Anne Soft Pack
- 02 50 21 American Red Cross Jacket for AED Little Anne
- 94 50 50 Remote Control
- 94 50 55 Programming Kit
- 94 50 90 Laerdal LINK Training Pads for AED Resusci Anne and AED Little Anne

Laerdal® AED Trainer 2

- 94 50 32 American Red Cross Soft Pack for AED Trainer 2
- 07-10900 FR/AED Trainer Training Pads (1 Set)

HeartStart® AED Trainer

- M5093A Adult Replacement Pads (1 Set)
- M5094A Infant/Child Replacement Pads (1 Set)
- M5087A Carry Case

AED Resusci Anne®

- 15 35 00 Soft Pack/Training Mat (for Torso Models)
- 11 01 00 Hard Case (for Full Body)
- 04 20 00 Hard Case (for Torso)

Limb Modules

- 31 03 00 Arms/Soft Legs with Soft Pack
- 31 20 00 Rescue Module with Soft Pack
- 31 20 50 First Aid Trauma Module with Soft Pack
- 31 20 09 Soft Pack for Arm/Leg Modules
- 31 03 21 Fastening Bolt, Arms (Pkg. 2)
- 31 03 41 Fastening Bolt, Lower Body

Output Devices

- 31 50 00 01 SkillReporter™ Complete (can only be used with Resusci Anne SkillReporter)

Training Supplies

- 18 39 10 Training Mat
- 31 03 51 Jacket
- 31 03 52 Trousers
- 31 50 50 Paper for SkillReporter (5 Rolls)

For CPR breathing barriers, see pages 13-16.



Laerdal® Face Shield CPR Barrier

Laerdal's most value-priced, convenient-to-carry, single-use CPR barrier

- Plastic sheet and 3M™ Filtrete™ hydrophobic filter help to reduce risk of blood or saliva contact between rescuer and victim
- Bi-directional filter and imprinted graphic directions for use

- 46 00 06 Laerdal Face Shield (Pkg. 50)
- 46 99 06 Laerdal Face Shield Key Chain (Pkg. 25)
- 46 00 14 Laerdal Face Shield Refills (Pkg. 25)



Bi-directional
3M™ Filtrete™
hydrophobic filter



Laerdal® Paediatric Pocket Mask®


The most widely used CPR mask in the world is now designed with infants and children in mind

- One-way valve acts as an additional barrier separating rescuer and victim
- Easy assembly for quick application
- Can be used in CPR instruction courses to teach mouth-to-barrier ventilation
- Distance from victim allows visualization of chest rise and fall
- Soft-pack carry case with belt strap and clip
- All components are latex-free

- 82 00 50 Laerdal Paediatric Pocket Mask w/Gloves and Wipe in Blue/Yellow Soft Pack
- 82 00 52 Laerdal Paediatric Pocket Mask w/Gloves and Wipe in Red/Black Soft Pack



Note: Wipes and gloves are expiration dated for user protection.




Laerdal® Pocket Mask™

Its one-way valve and 3M™ Filtrete™ hydrophobic filter give effective ventilations and outstanding user protection. Others may look like the Laerdal Pocket Mask, but there's only one original!

The Pocket Mask CPR barrier combines a low-resistance one-way valve with a disposable 3M Filtrete hydrophobic filter to help prevent the passage of liquids and secretions. Valve and filter act as barriers separating rescuer and victim.

- Reusable mask is easy to clean
- Increases distance from victim
- Pre-inflated cuff for ease of application, quick and effective seal
- Oxygen inlet version facilitates delivery of supplemental oxygen and includes latex-free head strap
- Can be used on adult, child, and infant
- Replaceable, disposable one-way valve and filter
- Latex-free
- Complete with latex-free and powder-free nitrile gloves and antimicrobial wipe, for the best available protection
- Versatile and attractive packaging options



Both the One-Way Valve and the Disposable Filter can be used with all previous versions of Laerdal Pocket Mask devices.



Note: Wipes and gloves are expiration dated for user protection.

Standard Versions

- 82 00 12 Pocket Mask with Gloves and Antimicrobial Wipe in White Hard Case (Each)
- 82 00 15 Pocket Mask in Soft Pouch with Gloves and Antimicrobial Wipe (Each)
- 82 00 14 Pocket Mask in Resealable Poly Bag (Each)

Oxygen Versions

- 83 00 06 Pocket Mask with Oxygen Inlet and Head Strap with Gloves and Antimicrobial Wipe in White Hard Case (Each)
- 83 00 15 Pocket Mask with Oxygen Inlet in Soft Pouch with Gloves and Antimicrobial Wipe (Each)
- 83 00 14 Pocket Mask with Oxygen Inlet and Head Strap in Resealable Poly Bag (Each)

Parts and Update Pack

- A** 82 06 10 Filter
- B** 82 04 10 One-Way Valve
- 82 07 00 Update Pack (1 Filter, 1 Valve)
- 83 01 10 Headstrap, Adjustable (Pkg. 10)

Easy, versatile, and durable carry case.



Resusci® Junior CPR Manikin with Water Rescue Kit

For water rescue training

The Water Rescue Update Kit converts your new or existing Resusci Junior into a Water Rescue Training Manikin. Installation does not affect the sanitation system.

- Naturally occluded airway can be opened with head-tilt/chin-lift and/or jaw-thrust techniques
- Carotid pulse simulation
- Convenient disposable airways to reduce cleanup time
- Removable/reusable faces
- Uses the same faces/airways as the Little Junior™ CPR training manikin



This lifelike CPR training manikin is specially designed for child resuscitation and water-rescue techniques.



The Resusci Junior Manikin with Water Rescue Update Kit

- 18 00 1201 Resusci Junior with Water Rescue Kit
- 18 38 20 Soft Pack Carry Case/Training Mat

For more information on Resusci Junior, see page 20.



- To update your existing Resusci Junior, we offer the following:
- 18 40 10 Water Rescue Update Kit
(Includes Leg Weights, Flotation Bladder, and Installation Instructions)

Resusci Anne® Basic and SkillGuide

Setting the standard for realistic CPR training

The Resusci Anne Basic and Resusci Anne SkillGuide manikins offer you an excellent starting point for delivering quality CPR training. The addition of the SkillGuide Performance Indicator provides immediate and objective feedback on key CPR elements.

- Anatomically correct landmarks and realistic compressions/ventilations for learning correct CPR techniques
- Correct head-tilt/chin-lift and/or jaw-thrust techniques must be performed for ventilation.



Resusci Anne CPR Manikin Includes:

Full-Body Manikin, 4 Decorated Resusci Manikin Faces, 3 Disposable Manikin Airways, Hard Carry Case or Soft Pack Carry Case, Track Suit (Fu I-Body Models) or Jacket (Torso Models), Directions for Use, and 50 Manikin Wipes.

Resusci Anne full-body comes standard with non-articulating arms and legs. First aid trauma modules available, see page 23.



Models listed as "SkillGuide" include the SkillGuide unit.


- 31 00 05 01 Resusci Anne Torso Basic w/Hard Case
- 31 00 06 01 Resusci Anne Torso Basic w/Soft Pack Carry Case
- 31 00 35 01 Resusci Anne Full-Body Basic w/Hard Case
- 31 00 15 01 Resusci Anne Torso SkillGuide w/Hard Case
- 31 00 16 01 Resusci Anne Torso SkillGuide w/Soft Pack Carry Case
- 31 00 45 01 Resusci Anne Full-Body SkillGuide w/Hard Case

Resusci Anne® SkillReporter™

Ideal for professional rescuer CPR training

The SkillReporter Resusci Anne manikin allows instructors to provide high-level CPR training for professional rescuers where optimal results are required. The intuitive and compact SkillReporter module has a built-in printer for documenting student performance against standard CPR protocols.

- Students get real-time feedback via LED display and/or graphical printed output
- Provides objective statistical report on accuracy in meeting standard CPR intervention protocols
- Has capability of logging CPR intervention sequence without the SkillReporter unit being connected, a very useful feature in scenario-based training



It is possible to perform and record the entire CPR sequence without the SkillReporter attached by using the manikin's Power and Log switches. Connect the SkillReporter unit later to download and report student's performance.



Resusci Anne CPR Manikin Includes:

Full-Body Manikin, 4 Decorated Resusci Manikin Faces, 3 Disposable Manikin Airways, Hard Carrying Case or Soft Pack Carry Case, Track Suit (Full-Body Models) or Jacket (Torso Models), Directions for Use, and 50 Manikin Wipes. Models listed as "SkillReporter" include the SkillReporter unit.

Optional articulating rescue and first aid modules can be added to perform realistic full-body scenario-based training. See page 23.



- 31 00 2501 Resusci Anne Torso SkillReporter
- 31 00 5501 Resusci Anne Full Body SkillReporter

Resusci® Baby

Quality CPR training with real-time feedback

The most realistic and convenient infant CPR training manikin available.

- Head-tilt/chin-lift is required to open the airway
- Brachial pulse simulation
- Economical, disposable airways for quick and easy cleanup
- Removable, reusable faces
- Optional SkillGuide™ CPR Performance Indicator provides immediate feedback on ventilation duration and volume, compression depth, and finger placement



Resusci Baby Manikin Includes:

Infant Manikin, Hard Plastic Carry Case, 3 Baby Faces, 6 Baby Airways, 6 Resusci Manikin Wipes, Pulse Bulb, Tubing, and Directions for Use.

- 14 00 1001 Resusci Baby Manikin
- 14 00 1101 Resusci Baby Manikin with SkillGuide Indicator

Model listed as 'SkillGuide' includes the SkillGuide unit.

Resusci® Junior

A child member of the Resusci family

This lifelike CPR training manikin is specially designed for child resuscitation and rescue techniques.

- Naturally occluded airway can be opened with head-tilt/chin-lift and jaw-thrust techniques
- Carotid pulse simulation
- Convenient disposable airways to reduce cleanup time
- Removable/reusable faces
- Uses the same faces/airways as the Little Junior™ CPR Training manikin
- Optional SkillGuide™ CPR Performance Indicator provides immediate feedback on ventilation duration and volume, compression depth, and finger placement



Resusci Junior Manikin Includes:

Pediatric manikin, Hard Plastic Carry Case (except where noted), Training Mat, Track Suit, 4 Junior Faces, 4 Junior Airways, 12 Resusci manikin Wipes, and Directions for Use.

Model listed as 'SkillGuide' includes the SkillGuide unit.

Model listed as Water Rescue includes the Water Rescue Update Kit accessory.

- 18 00 10 01 Resusci Junior Manikin with Training Mat
- 18 00 13 01 Resusci Junior Manikin with Training Mat
- 18 00 11 01 Resusci Junior Manikin with SkillGuide and Training Mat

Resusci Junior can also be used for water-rescue training. See page 17.



Laerdal® Airway Management Trainer

Lifelike airway develops skill proficiency

Realistic practice is the key to developing proficiency in airway management skills. The Laerdal Airway Management Trainer's lifelike upper torso and head simulate real-world complications when practicing a variety of intubation, ventilation, and suction techniques.

- Realistic anatomical features allow demonstration of Sellick Maneuver and laryngospasm
- Improved airway allows insertion of Laryngeal Mask Airway and Combitube®
- Practical training in clearing an obstructed airway and suctioning of liquid foreign matter
- Provides visual inspection of lung expansion and auscultation of breath sounds
- Simulates stomach inflation and vomiting situation
- Airway Demonstration Model is standard with each trainer



25 00 0033 Laerdal Airway Management Trainer


Laerdal Airway Management Trainer Includes:

Manikin on Sturdy Board, Airway Demonstration Model, Cleaning Kit, Airway Lubricant, Carry Case, and Directions for Use.




Suction

Unique possibilities for practicing suction from the oral cavity, oropharynx or nasopharynx, orotrachea or nasotrachea via endotracheal tube.



Ventilation

Permits realistic manual maneuvers to open the airway. Suitable for practicing bag-valve-tube/mask ventilation.



Intubation

Allows practicing tracheal, pharyngeal, and esophageal intubation. Oral or nasal.



Improved airway allows insertion of Laryngeal Mask Airway and Combitube®

AED Resusci Anne® SkillReporter™ and SkillGuide™

The standard in realistic CPR/AED training

The AED Resusci Anne combined with our Laerdal AED Trainer 2 allow the instructor to control the training scenario and the students to train to perfection. The system requires proper pad placement before the scenarios continue.

- LEDs in the manikin's chest confirm ideal pad placement
- Added instructor flexibility with the optional Laerdal AED Trainer 2 remote control allows selecting, pausing, resuming, and overriding the scenario
- Events logged in this system include:

- ▶ Power on
- ▶ Pads attached in proper place
- ▶ Loose electrodes
- ▶ No shock indicated
- ▶ Motion artifact
- ▶ Low battery
- ▶ AED error
- ▶ Pause scenario



- 32 00 6001 AED Resusci Anne Torso SkillGuide
- 32 00 7001 AED Resusci Anne Torso SkillReporter
- 32 00 8001 AED Resusci Anne SkillGuide
- 32 00 9001 AED Resusci Anne SkillReporter

All Resusci Anne manikins Include:
 4 Decorated Resusci Manikin Faces, 3 Disposable Manikin Airways, Hard Carrying Case or Soft Pack Carry Case, Track Suit (Full-Body Models) or Jacket (Torso Models), and Directions for Use. Models listed as "SkillGuide" include the SkillGuide unit. Models listed as "SkillReporter" include SkillReporter unit with 1 roll of printer paper.

Resusci Anne® First Aid/Trauma Manikin Accessories

- Weighted, articulated arms and legs with simulated wounds can easily be attached to the manikin within a few seconds. No tools are needed
- The arms and legs feature a selection of simulated burns, cuts, and fractures
- Allows you to train according to Lifesaving First Aid or Basic Trauma Life Support protocols
- The wounds may alternatively serve as distracting elements for added realism in CPR and CPR-Defibrillation training
- Limb articulation allows for realistic victim handling and immobilization training. New short hairstyle allows for realistic application of extrication collar
- Penetrating bleeding wounds with fractured femur (bleeding control leg) allow you to train control of bleeding
- Optional IV arm available



- 09 00 50 Bleeding Control Leg
- 31 20 25 IV Arm
- 31 20 50 First Aid/Trauma Module with Soft Pack
- 31 20 00 Rescue Module with Soft Pack

Bleeding Control Leg 09 00 50



Penetrating bleeding wounds with fractured femur allow training of bleeding control and limb stabilization.

IV Arm 31 20 25



IV arm incorporates IV skills into your training scenarios.

First Aid/Trauma Modules 31 20 50



Arms and legs feature a selection of simulated burns, cuts, and fractures.

Rescue Module 31 20 00



Weighted articulated arms and legs allow for training of rescue and extrication techniques combined with CPR or CPR-Defibrillation.

Resusci Anne Continued

Airway and Faces

Resusci Anne®

- 15 22 50 Airways 700-1000 ml (Pkg. 24)
- 15 26 00 Airways 400-600 ml (Pkg. 24)
- 31 02 10 Manikin Faces (Pkg. 6)
- 31 02 20 Manikin Faces Decorated (Pkg. 6)

Resusci® Junior

- 18 30 10 Junior Faces (Pkg. 6)
- 18 32 10 Junior Airways (Pkg. 25)
- 18 32 11 Junior Airways (Pkg. 100)

Resusci® Baby

- 14 37 00 Resusci Baby Airways (Pkg. 24)
- 14 37 01 Resusci Baby Airways (Pkg. 96)
- 14 36 00 Resusci Baby Faces (Pkg. 6)

Carrying Cases

Resusci Anne

- 15 35 00 Soft Pack/Training Mat (for Torso Models)
- 11 01 00 Hard Case (for Full-Body)
- 04 20 00 Hard Case (for Torso)
- 31 20 09 Soft Pack for Arm/Leg Modules

Resusci Baby

- 14 40 10 Hard Case

Miscellaneous

Airway Management

- 25 20 90 Airway Lubricant
- 15 24 00 Manikin Wipes (Pkg. 50)
- 15 24 01 Manikin Wipes (Pkg. 1,200)
- 15 12 0103 Manikin Face Shields (6 Rolls of 36 Each)

Resusci Anne®

Output Devices

- 31 50 00 01 SkillReporter™ Complete (can only be used with Resusci Anne SkillReporter)

Limb Modules

- 31 03 00 Arms/Soft Legs with Soft Pack
- 31 20 00 Rescue Module with Soft Pack
- 31 20 50 First Aid Trauma Module with Soft Pack
- 31 03 21 Fastening Bolt, Arms (Pkg. 2)
- 31 03 41 Fastening Bolt, Lower Body

Training Supplies

- 18 39 10 Training Mat
- 31 03 51 Jacket
- 31 03 52 Trousers
- 31 50 50 Paper for SkillReporter (5 Rolls)

Resusci Anne with SkillGuide™ Option

- 15 40 00 Update Kit, SkillGuide

Resusci Anne with Defibrillation Option

- 32 10 00 Defibrillation Update Kit, Basic/SkillGuide
- 32 10 10 Defibrillation Update Kit, SkillReporter

Resusci Anne with Laerdal Link Technology Option

- 32 10 60 AED w/Laerdal Link Technology Update Kit, SkillGuide
- 32 10 70 AED w/Laerdal Link Technology Update Kit, SkillReporter

Resusci® Junior

Update Kit for Older-Style

Resusci Junior Manikins

- 18 40 10 Resusci Junior Water Rescue Update Kit

Converts your new or existing Resusci Junior (Catalog numbers: 18001301 and 18001001) into a manikin that can be used for water-rescue training. Installation does not affect the sanitation system. Includes leg weights, a bladder for flotation, and installation instructions.

Airway Bags

- Designed for older-style Resusci Junior airway systems without disposable airways.
- 18 32 04 Resusci Airway Bags (10 Rolls of 50 Each)

Laerdal® PC SkillReporting System

Quality CPR training and evaluation tool

The Laerdal PC SkillReporting System offers an approach to CPR or CPR-D training with Resusci Anne® by integrating the benefits of a desktop or laptop PC into your training. The system includes an extensive student/course administration functionality. When accurate documentation of student performance is a requirement or you would like to run and control multiple manikins from one PC, this is the perfect tool for you.



317000 Laerdal PC SkillReporting System

(includes: Program CD and Interface Cable. PC not included.)

With the Laerdal PC SkillReporting System you can operate multiple Resusci Anne SkillReporter manikins with one PC.

One of the following Resusci Anne® models is required:

Resusci Anne SkillReporter



Resusci Anne Advanced SkillTrainer



Resusci Anne Simulator



Thomas™ Tube Holder

State-of-the-art improvement over tape-fastening methods

- Accommodates multiple size ET tubes as well as Combitubes and LMAs
- Bite block protects victim's teeth and gums while keeping the ET tube open
- BVM does not need to be removed to attach the Thomas Tube Holder
- Allows for in-use suctioning of the oropharynx
- Disposable – one-time use



600-10000 Thomas Adult Tube Holder (only sold in cartons of 25)

600-20000 Thomas Pediatric Tube Holder (only sold in cartons of 25)

Trauma Make-Up Kit

- Create "shock," burns, and other wounds
- Simulated blood concentrate
- Six colors of grease paint
- Modeling putty wax
- Knox® Gelatin
- Skin adhesive
- Mist brushes
- Atomizer
- Baby powder
- Tongue depressors
- Carrying case (not pictured)



Trauma Make-Up Kit Includes:

Simulated Blood Concentrate (4), Grease Paint - Blood Red (2), Grease Paint - Red, Grease Paint - Blue, Grease Paint - White, Grease Paint - Ebony, Modeling and Putty, Wax (2), Knox Gelatin, Skin Adhesive, Artist Brushes (2), Atomizer, Baby Powder, Tongue Depressors (5), Carry Case

275-10001 Trauma Make-Up Kit

Laerdal® Silicone Resuscitator

The reusable resuscitator with design and performance based on over 35 years of experience

- New masks provide a better seal and are easier to clean
- Multi-function mask cover and new O₂ intake valve
- See-through mask design enables clear view of victim at all times
- Resuscitator is easy to use. Offers the sensitivity that signals changes in lung compliance
- Choice of sizes
- Superior bag re-expansion, resistant to high/low temperatures
- Easy cleaning, autoclavable to 276°F
- Latex-free



- 87 00 5233 Adult Standard, w/Mask, in Carton
- 86 00 5233 Pediatric Standard, w/Mask, in Carton
- 85 00 5133 Preterm Complete, in Carton

The BAG™ II Disposable Resuscitator

Performance and value... it's in The BAG

- Adjustable hook-and-loop handle for all hand sizes and maximum comfort
- The texture-ribbed body facilitates an easier grip, minimizes fatigue, and provides an excellent feel of victim lung compliance
- Shaped to maximize ventilation volume
- Performs well in extreme temperatures: 0° to 122°F
- Disposable Face Mask features an ultra-soft airtight cushion seal
- A 10 ft. connecting tube provides easy access to oxygen source
- Optional Disposable PEEP Valve is adjustable and position independent




- 84501 I The BAG Adult Disposable BVM (individual)
- 84502 I The BAG Child Disposable BVM (individual)
- 84503 I The BAG Infant Disposable BVM (individual)

V-VAC™ Manual Suction Unit


The all-in-one manual suction system

Get a good grip on handling even the toughest suction scenario with Laerdal's V-VAC Manual Suction Unit. It's ready to go, easy to use, and up to the task. Ideal for First Responder emergency kits and as back-up suction, it meets all your suction needs, infant through adult. One-handed operation, one-piece disposable cartridge, and unique no-clog suction tip make the V-VAC the manual suction unit of choice!

- Adjustment for high- and low-pressure suctioning
- Efficient: Maximum suction exceeds 380 mmHg. Peak air flow over 70 liters per minute
- Features a unique no-clog suction tip (intake valve)
- Versatile catheter connection
- Disposable one-piece cartridge provides a clean, new pump, valve, and suction tube
- Reusable handle is immersible for cleaning



V-VAC unit fits into many emergency kits. A separate carrying bag is also available.



The soft tip of the intake valve has deep cutouts to prevent tube occlusion, and also soft-tissue injury, when the tip is held against the back of the throat. The tri-cut valve resists back flow of suctioned matter.

- 98 50 00 V-VAC Starter Kit
- 98 51 00 V-VAC Training Kit



V-VAC Training Kit for Practical Use Includes:
Handle, 1 Reusable Training Cartridge w/Plug for Emptying and Cleaning, Catheter, and Connectors.

V-VAC Starter Kit (shown) for Patient Use Includes:
Handle, 2 Replacement Cartridges, Short Suction Catheter w/Adapter Tip, Double Male Connector, and Directions for Use.

Accessories

- | | | | |
|----------|-----------------------|----------|----------------------------------|
| 98 50 01 | Replacement Cartridge | 98 50 02 | Adapter Tips (Pkg. 4) |
| 98 51 01 | Training Cartridge | 98 50 03 | Double Male Connectors (Pkg. 10) |
| 98 52 00 | Carrying Bag | 98 50 04 | Catheters (8" / 18 F, Pkg. 4) |

Stifneck® Select™ Adjustable Collar

Savvy providers select quality with the Select. It's simple symmetrical sizing plus safety!

The Stifneck Select collar simplifies your job, saves space, and provides outstanding immobilization. The quickest and easiest-to-use adjustable collar on the market, all you have to do is measure the patient, select the size, adjust, and lock!

By packing the Select Collar, you'll have more room in your bag for other equipment.

This collar fits all adult patients and delivers the same superior immobilization you've come to expect from Stifneck collars.

- Reduces inventory costs
- Uses the same sizing and application method of the original Stifneck collars
- Adjustment tracks ensure symmetrical alignment of size
- Easy access for pulse checks, advanced airway procedures, and visualization through oversized trachea hole
- Room for large fingers to slide through rear panel opening for palpation



- 98 00 10 Stifneck Select Collar
- 98 19 20 Stifneck Select Training Video (NTSC)

Stifneck® Pedi-Select™ Extrication Collar

- Superior immobilization to fit your smaller patients
- Pre-assembled chin piece ensures stability
- Optimal visualization for assessment, pulse checks via the front and back openings

- 98 00 20 Stifneck Pedi-Select Collar



Stifneck® Extrication Collars*


The classic collar against which all others are measured

Refined for over 20 years, the original Stifneck extrication collar remains the standard of care for suspected cervical spine injuries. Though there are other collars available, the Stifneck is the only collar that measures up.

- Outstanding motion restriction without hyperextension or limiting access to airway
- Provides optimal size selection with simple, accurate finger-sizing method and color codes
- Fits entire patient population with four adult and two pediatric sizes
- Independent, peer-reviewed studies support safety and efficacy
- Radiolucent, MRI and CT scan compatible
- Latex-free



- A** 98 06 00 Tall Collar (Adult)
- B** 98 05 00 Regular Collar (Adult)
- C** 98 04 00 Short Collar (Adult)
- D** 98 03 00 No-Neck Collar (Adult)
- E** 98 02 00 Pediatric Collar
- F** 98 01 00 Baby No-Neck™ Collar



Pediatric and Baby No-Neck collars offer a rear panel opening for palpation, visualization, and ventilation.

SpeedBlocks® Head Immobilizer


The first high-performance reusable at the cost of a disposable

Designed for durability and performance, but priced for disposability, SpeedBlocks really give you the best of both worlds.

SpeedBlocks offers a simple, four-step process for application, a unique, multi-axis blocks system to accommodate all head shapes and injury positions, and direct access to patient's ears for better communication and inspection.

The low-cost blocks with disposable foam pads and head straps are so effective that the inconvenience of taping the patient can be eliminated.

- Disposable head straps and block pads allow for easy and low-cost replacement
- Reliable handle-locking mechanism with quick-release locks ensures confident application
- Blocks provide contoured fit for patients aged 2 and up
- Universal base and blocks with replaceable foam pads for easy, worry-free cleaning
- Blocks slide for easy multi-axis adjustment
- Blocks weigh only 5.2 oz each
- Radiolucent, MRI and CT scan compatible
- Made of HDPE durable plastic and non-absorbing padded blocks
- Latex-free



- 98 30 90 SpeedBlocks Starter Pack
(Includes 1 Block Set, Universal Base, and 1 set of Pads)
- 98 30 92 SpeedBlocks Blocks Set
- 98 30 60 SpeedBlocks Universal Base
- 98 30 80 SpeedBlocks Training Video
- 98 30 96 SpeedBlocks Strap and Pad Replacement Set
(Includes 5 sets of Block Pads, Head Straps, and Base Pads)

Sta-Blok™ Head Immobilizer

Our disposable, all-weather foam HID! Easy access, easy adjustment, easy on the budget. Foam headrest helps elevate spine to neutral position while providing victim comfort

- Easy to re-adjust with Velcro® padded head strap
- Includes "Severe Weather/Helmet" strap to provide extra stability
- Can fit virtually any patient on any spineboard
- Ears and head accessible without removing the device for ongoing patient assessment



700-0001 Sta-Blok Head Immobilizer
(Pkg. 30)

HeadBed™ Head Immobilizer

All the strength of corrugated board in all-weather secured, disposable immobilizer

- Disposable, so no equipment recovery costs or cleanup
- Effective in various weather conditions
- Outstanding patient access for effective assessment without compromising immobilization
- Fits adults, children and infants
- Fits most spine boards
- Radiolucent, MRI and CT scan compatible
- Our patented Dual Adhesive Strap sticks aggressively where needed and releases where necessary
- Latex-free



98 20 00 HeadBed II HID (each)

PadPack™ Alignment Pads

- Assist in the maintenance of neutral alignment:
Adult – Occipital padding, Child – Torso padding
- Simple and easy to use



98 21 00 PadPack Alignment Pads (Pkg. 10)

PadPack Pad

BaXstrap™ Spineboard

The board that meets your most demanding needs

As a rescuer, you want your spineboard to be comfortable for your victim, rigid, and lightweight. The BaXstrap spineboard provides this and is one of the lightest boards on the market. The contoured design allows the board to nest compactly. This unique feature also offers patients more comfort, making it less likely that they would move around and aggravate their injuries.

- Engineered for maximum strength and designed for both patient and rescuer comfort
- No seams make it easy to clean and decontaminate
- Pins are molded into the board, maximizing their strength
- Includes strap holes for pediatric patients, eliminating the need to use towels to fill the extra space between patient and strap
- Can be used in water rescue
- Durable, molded-in custom labeling available

98 25 00 BaXstrap Spineboard

98 25 12 Premium Strap with sewn loop ends (each)

98 25 13 Premium Strap with speed clip ends (each)



Lifting is made easier for the rescuer with raised handholds designed to maximize ground clearance.

Can be easily stored. Three boards stacked together are only 5 1/2 inches thick.

Laerdal® Pedi-Pad Spineboard Pad

It's the foam, anti-skid kid disposable pad that provides proper alignment in any spineboard! Pedi-Pad works for kids

- Places the cervical spine in neutral alignment by elevating the pediatric patient's shoulders
- Economical method to stabilize children 8 years old and younger

701-00001 Laerdal Pedi-Pad Spineboard Pad (Pkg. 6)



It is recommended that you fill out the enclosed order form first (see page 35), indicating the catalog numbers, quantities, product descriptions, and prices. Once the order form is complete, please contact your participating American Red Cross chapter (see address below) to place your order.

Your chapter's/station's address, telephone, and fax numbers:

Helpful Hints

1. When ordering, please note the package quantities of various parts and accessories.
2. Call your participating American Red Cross unit if you have any questions.

Payment Methods

1. When ordering from the American Red Cross, make checks or money orders payable to your participating American Red Cross unit.
2. Contact the participating American Red Cross unit for information on credit card orders.



**American
Red Cross**



Laerdal®

helping save lives

HeartStart® is a registered trademark of Philips Medical Systems. 3M™ and Filtrete™ are trademarks of 3M. Velcro® is a registered trademark of Velcro Company. Combitube® is a registered trademark of Kendall, a division of Tyco Healthcare. Knox® is a registered trademark of the Knox Company.

©2009 Laerdal Medical Corporation. All rights reserved. Printed in the USA. #4697
PRO-ML02-0166.A