

Clinical skills training for hospital professionals

Now with  
ECGs!



## Nursing Anne

*The manikin for the development of clinical skills of hospital professionals*

Nursing Anne is designed for scenario based training for the care and management of a wide variety of hospital patients. Nursing Anne is an efficient, effective, flexible manikin for clinical training in women's health, obstetrics, post-partum, general patient or wound assessment and care.

- Educationally effective for hospital training **targeting key skills of hospital care**, including women's health, obstetrics, post-partum care, wound assessment and care and general patient care
- Preprogrammed scenarios provide **standardized training** while customizable scenarios and real-time instructor control allows adaptation to meet individual student needs
- **Flexible manikin platform** allows multiple accessory modules to be added, including nursing wounds, breast exam, fundus assessment and mastectomy
- Head with anatomical landmarks, trachea and esophagus, along with simulated lungs and stomach, allow the practice of many procedures, including NG, OG, tracheal care and suctioning
- Normal and abnormal heart, breath and bowel sounds, and fetal heart tones for auscultation
- Instructor-controlled blood pressure arm allows for realistic palpation and auscultation. Systolic and diastolic pressures, auscultatory gap and volume are variable



# Nursing Anne Product Specifications

## Patient Care

Irrigation of the eye and ear (simulated)  
 Nasal packing  
 Mouth and denture care procedures  
 Insertion and suctioning of oropharyngeal and nasopharyngeal airways  
 Tracheostomy care and tracheal suctioning  
 Various oxygen delivery procedures  
 Manually generated carotid pulse  
 NG tube insertion, care and removal  
 Gastric lavage and gavage  
 Subclavian IV, various ostomy and chest tube openings (placement only)  
 Injection sites include bilateral deltoid, bilateral thigh, gluteal and ventrogluteal  
 Articulating IV arm with replaceable skin and infusible vein system allows peripheral intravenous therapy and site care  
 - Venipuncture possible in the antecubital fossa and dorsum of the hand  
 - Accessible veins include median, basilic and cephalic  
 Interchangeable male and female genitalia with connector valves for urinary and colon reservoirs  
 - Complete urinary catheterization  
 - Will retain indwelling or straight catheter  
 - Enema procedures may be performed using fluid for realistic return  
 Full range of motion for realistic patient handling

## Blood Pressure / Pulses – VitalSim™

Articulating blood pressure arm for auscultated and palpated blood pressure simulation  
 Korotkoff sounds synchronized with programmable ECG  
 Korotkoff sounds volume control  
 Systolic and diastolic pressure may be set individually in steps of 2 mmHg  
 Systolic 0–300 mmHg, Diastolic 0–300 mmHg  
 Auscultatory Gap, on/off feature  
 Pressure accuracy +/- 2 mmHg  
 Calibrate function to adjust pressure sensor and cuff gauge  
 Palpated pulse simulator controls brachial and radial pulses  
 - Pulses synchronized with programmable ECG  
 - Pulse strengths dependent on BP

325-05050 Nursing Anne (VitalSim™ Capable)  
 200-10001 VitalSim™ Control Unit

**Nursing Anne VitalSim™ Capable Includes:**  
 Full-body Adult Manikin, Hospital Gown, Manikin Lubricant, Assembly Tool Kit and Directions for Use

## Sounds – VitalSim™

Heart sounds synchronized with programmable ECG  
 Lung sounds synchronized with breathing rate, 0–60 BPM.  
 Individual lung or bilateral sound selection  
 Normal and abnormal bowel sounds  
 Pre-recorded or live feed voice sounds

## Logging / Scenario Functions – VitalSim™

The PC editor allows for creating and editing of scenarios  
 Up to 10 scenarios and 25 logs held in the VitalSim Control Unit  
 Logs may be downloaded to PC for after-action review

## Cardiac Capabilities **New!**

Now with ECGs  
 1400+ rhythm variations for monitoring and interpretation

## Optional Accessories



325-00750  
 Breast Exam Module



325-00650  
 Mastectomy Module



325-00450  
 Fundus Skills and Assessment Module



325-00550  
 Wound Care and Assessment Set



For more product information, please visit [www.laerdal.com](http://www.laerdal.com)

**For more information, please call (800) 431-1055 (US), (888) 523-7325 (Canada), (800) 567-9987 (Quebec)**

**USA and Latin America:**  
 Laerdal Texas  
 P.O. Box 38-226 FM 1116  
 Gatesville, Texas 76528  
 (800) 433-5539 or  
 (254) 865-7221  
 FAX (800) 865-8011  
[www.laerdal.com](http://www.laerdal.com)

Laerdal Medical Corporation  
 167 Myers Corners Road  
 Wappingers Falls, NY 12590  
 (800) 431-1055 or  
 (845) 297-7770  
 FAX (800) 227-1143  
[www.laerdal.com](http://www.laerdal.com)

**Canada:**  
 Laerdal Medical Canada Ltd.  
 151 Nashdene Rd., Unit #45  
 Toronto, ON M1V 4C3, Canada  
 (888) LAERDAL  
 (416) 298-9600  
 FAX (416)-298-8016  
[savelives@laerdal.ca](mailto:savelives@laerdal.ca)

