

Material Information Sheet

Last updated: 19.07.2007

Catalogue number	Product name
980600	Stifneck Tall Collar (Adult)
980500	Stifneck Regular Collar (Adult)
980400	Stifneck Short Collar (Adult)
980300	Stifneck No-Neck Collar (Adult)

Part/ Group of parts	Material	CAS-no (if available)	Comments
Neckband	High-density polyethylene		Yellow Color
Chin Support	High-density polyethylene		Blue Color
Rivet, Male	Nylon		
Hook Collardot	Nylon		
Hook and Loop	Nylon		
Foam pads	Low-Density Polyethylene	9002-88-4	Biocompatible according to ISO 10993-1 "Safety for use with humans, skin contact <24 h"
Carry Bag Stifneck	Polyester fabric with PVC backing		

This product complies with:
California Proposition 65 (20.04.2007)

This product contains no latex.

Author:

Eivind Tønnessen
Technical Product Specialist R&D
Laerdal Medical