

Resusci Anne Q-CPR/D training system


Quality CPR combined with the use of HeartStart MRx (manual or automatic) for maximum realism

Resusci Anne Q-CPR/D training system allows for training with the Q-CPR function on the HeartStart MRx defibrillator. The system allows for defibrillation, Q-CPR and airway management training. The airway head is designed for controlled ventilations through use of endotracheal tube, laryngeal mask or bag-valve-mask ventilations.

The Q-CPR Training Module provides ECG to the defibrillator and accepts shocks. It also provides ventilation data to the HeartStart MRx.

ECG and shock response is controlled from the Micro-HeartSim.

With the use of the optional SkillReporter printer or the PC SkillReporter comprehensive debriefing documentation can be generated.



199-10050 Q-CPR Training Module. *Must be bought separately.*

The Q-CPR Training Module can be connected to any Resusci Anne CPR/D training manikin and to the special Resusci Anne Q-CPR / D manikin.

The module connects directly to the HeartStart MRx defibrillator contact.

199-00050 Resusci Anne CPR-D Torso with Airway Head.

199-05050 Resusci Anne CPR-D Torso with Airway Head.

Includes MicroHeartSim ECG simulator.

Consumables

25 20 90
Airway
Lubricant



315050
Paper for
SkillReporter
(pkg. 5)

