

Material Information Sheet

Last updated: 20.06.2007

Catalogue number	Product name
980012	Stifneck Select Collar, Green

Part/ Group of parts	Material	CAS-no (if available)	Comments
Chin support	High-Density Polyethylene		Green color
Neckband	High-Density Polyethylene		Green color
Rivet, Male	Nylon		
Loop fastener	Nylon		Black color
Foam pads	Low-Density Polyethylene	9002-88-4	Black color Biocompatible according to ISO 10993-1 "Safety for use with humans, skin contact <24 h"

This product complies with:
California Proposition 65 (20 04 2007)

This product contains no latex or phthalates

Author:

Eivind Tønnessen
 Technical Product Specialist R&D
 Laerdal Medical