

**The first high-performance reusable
at the cost of a disposable.**



SpeedBlocks® Head Immobilizer

Designed for durability and performance, but priced for disposability, the SpeedBlocks Head Immobilizer really gives you the best of both worlds.

The low-cost blocks with disposable foam pads and head straps are so effective that the inconvenience of taping the patient can be eliminated.

SpeedBlocks offer a simple four-step process for application, a unique multi-axis blocks system to accommodate all head shapes and injury positions, and direct access to patient's ears for better communication and inspection.

Low-cost, disposable...

- **Block & Base Foam Pads**
- **Velcro® Head Straps**

Superior patient access

Contoured fit for increased patient comfort

Multi-axis sliding adjustment

Reusable base and blocks



Laerdal™

helping save lives®

Product specifications

FEATURES	BENEFITS
Snap-in blocks	<ul style="list-style-type: none"> • Quick and easy application
Blocks slide for adjustment	
Wide lateral adjustment	<ul style="list-style-type: none"> • Patient stability • Allows re-adjustment with little risk of moving head • Can splint head as found • Fits all head shapes and positions
Head straps attached at top of blocks	
Contoured padded blocks	<ul style="list-style-type: none"> • Fits most helmets (use with tape) • Quick access • Can be instantly applied and removed • Easy to re-adjust • Not affected by rain or temperature
Low profile open design	
Large, low profile ear holes	<ul style="list-style-type: none"> • Fit device to patient (fits patients 2 years and up) • Patient comfort • Clear and easy patient access • Less intrusive around head, neck, face and ears
Velcro®-free cleanable base with disposable foam base pad	
Universal base	<ul style="list-style-type: none"> • Patient can hear more clearly • Better access and visualization of ears • Easily cleaned and disinfected, lowering cost • Can be used with most spineboards • Helps maintain neutral alignment
Low profile base	
X-ray translucent	<ul style="list-style-type: none"> • Stores attached to spineboard and ready for use • Patient remains immobilized until cleared or confirmed • Helps provide clear diagnosis
CT & MRI scan compatible	

With an even greater focus on Spinal Cord Injuries (SCI) and the possibility of using our “training programs”, you as an emergency care provider, can help reduce the chances of SCI. The Stifneck® Select™, SpeedBlocks and BaXstrap™ Spineboard are all designed with this in mind!



Item	Single	Case
SpeedBlocks® Starter Pack <i>(Includes 1 block set, 1 universal base and 2 block pads, 2 head straps, and 1 base pad)</i>	98 30 90	98 30 91 (Qty: 10 of 98 30 90)
SpeedBlocks® Blocks Set <i>(Includes 1 set of blocks)</i>	98 30 92	98 30 93 (Qty: 20 of 98 30 92)
SpeedBlocks® Universal Base <i>(Includes 1 base with base pad)</i>	98 30 60	98 30 70 (Qty: 10 of 98 30 60)
SpeedBlocks® Strap & Pad Replacement Set <i>(Includes 10 foam block pads, 10 head straps and 5 base pads)</i>	98 30 96	98 30 97 (Qty: 10 of 98 30 96)



For additional information or to order, contact your local Laerdal Authorized Distributor!

Laerdal Medical Corporation
167 Myers Corners Road, Wappingers Falls, NY 12590
(800) 431-1055 (845) 297-7770
FAX (800) 227-1143 www.Laerdal.com

Laerdal Medical Canada Ltd.,
151 Nashdene Rd., Unit #45 Toronto, ON M1V 4C3, CANADA
(416) 298-9600 (888) LAERDAL FAX (416) 298-8016
savelives@Laerdal.ca